

MEASURING TOXICITY

What you need:

- OP50
- NGM plates
- M9 buffer
- 15 mL tube
- 96-well microplate
- WMicroTracker device


Notes:

- Use M9 buffer supplemented with bacteria OP50 to avoid the static of micro plates.
- Use a saline buffer for tests no longer than 12 hours. The starved worms could affect the results.
- Perform at least three technical replicates and at least two biological replicates.


Protocol:

1. Grow synchronized populations of adult worms in seeding NGM plates (OP50).
2. Remove worms from plates using M9 buffer and transfer them in a sterile 15 ml tube.
3. Let the worms settle. Decant the supernatant taking care not to disturb the pellet.
4. Perform a wash with 5 ml of M9 buffer. Briefly shake or invert the tube.
5. Repeat the decantation step. Throw out the supernatant and add 3 ml of M9 buffer.
6. Count number of worms in 10 ul in triplicate and calculate the average.
7. Prepare a suspension to get [5 worms/10 ul]. Adjust volume in M9 buffer supplemented with bacteria OP50 (0.1 OD600 to final concentration).
8. Transfer 90 ul of worm solution to 96-well microplates using multichannel pipette. Let worms rest for 1 hour and measure basal activity for at least 30 minutes using WMicroTracker.
9. Add 10 ul of a 10x concentrated solution of chemicals to test. Include a control without compounds.
10. Register worm activity using wMicroTracker.

Data produced using this protocol:



Measuring KCl Concentration Effect. In this experiment we can observe the effect of increasing concentrations of KCl in young adults N2 *C.elegans*. The recording of 12h of activity shows an immediate effect of this salt on behavior and deleterious effect at doses higher than 3.75g/l.



Dose Response Of Aldicarb And Levamisole. In these experiments we can observe the kinetic and dose response to two Acetylcholine channels modulators: Aldicarb and Levamisole. In less than one hour a quantitative dose response effect is obtained.